



Q.VITEC

Erkennen, was möglich ist.

Q.VITEC Camera Manual

Q.VITEC GmbH
Hagenburger Str. 54,
31515 Wunstorf, Germany
Tel.: +49(0)5031-949 43 20
Fax: +49(0)5031-949 43 29
Web: www.qvitec.de



This manual and everything described in it are copyrighted. You may not copy this manual, in whole or part, without written consent of Q.VITEC GmbH (Q.VITEC).

Q.VITEC pursues a policy of continuous improvement of the design and performance of its products. Therefore we reserve the right to change the manual/product without notice. In no event will Q.VITEC be liable for direct, special, incidental, or consequential damage resulting from any defect in the product or its documentation, even if advised of the possibility of such damages.

Please direct support matters and technical questions to Q.VITEC.

© This description may not be reprinted or copied, in full or in part, without the express permission of Q.VITEC.

Q.VITEC GmbH
Hagenburger Str. 54,
31515 Wunstorf, Germany

Important symbols

One or more of the following symbols may be used in this documentation:



Warning!

The warning triangle indicates especially important safety instructions. If they are not adhered to, the results could be fatal or critical injury.



CAUTION

Indicates that you should proceed with caution. Failure to do so may result in injury or significant damage to instruments or their contents, e.g. data.



NOTE

Contains important additional information.



EXAMPLE

Contains an illustrative example of the previous text section.



PROCEDURE

Indicates that a step-by-step procedure follows.



REFERENCE

Indicates where you can find additional information on the subject at hand.



◆ KEY POINTS

Summarizes key points in a concise manner.



◆ SHORTCUTS

Provides helpful keyboard shortcuts.



◆ EXPLANATION

Provides brief explanation of a function, e.g. why or when you should use it.

➡ next page

Indicates that the text will be continued on the next page.

The manual uses the following conventions to indicate elements from the user interface or the keyboard:

| | |
|---------------------|--|
| "Data field" | Data field entries and option names are rendered in quotation marks. |
| [Button] | Buttons are indicated by square brackets. |
| <Key> | Keys are indicated by pointed brackets. |

Table of contents

- 1. Camera Specifications 7**
 - 1.1 Introduction8
 - 1.1.1 Camera types.....8
 - 1.1.2 Connecting Cameras to the Imagechecker11
 - 1.2 2-Mega-Pixel Digital Camera ANPVC 1210 / 2260..... 12
 - 1.2.1 Time Diagram for Image Grabbing12
 - 1.2.2 Dimensions14
 - 1.2.3 Specifications14
 - 1.2.4 Spectral Response15
 - 1.3 Quad-Speed Digital Camera ANPVC 1040 16
 - 1.3.1 Time Diagram for Image Grabbing16
 - 1.3.2 Dimensions18
 - 1.3.3 Specifications18
 - 1.3.4 Spectral Response19
 - 1.4 5-Mega-Pixel Digital Camera ANPVC 1510 20
 - 1.4.1 Time Diagram for Image Grabbing20
 - 1.4.2 Dimensions22
 - 1.4.3 Specifications23
 - 1.4.4 Spectral Response23
 - 1.5 Partial Scan Mode24
- 2. Index..... 27**
- 3. Record of Changes 28**

Chapter 1

Camera Specifications

1.1 Introduction

This chapter deals with the four different standard Camera Link types of cameras available for the Imagechecker. You will find information on how to connect a camera to the Imagechecker as well as a detailed description of each camera type.

You will find the following information:

- available camera modes (progressive scan, partial scan)
- how to use a flash
- camera hardware data (dimensions)
- tables with specifications.

Please read the information on the different camera types and their configuration carefully. The right camera choice optimizes results in processing speed and image resolution.

1.1.1 Camera types

The camera types listed below are available as standard types. Select the type most suitable for your application.

| | 2-Mega-pixel Gray digital camera (ANPVC 1210D) | 2-Mega-pixel Gray digital camera (ANPVC 1211D) → For details specification contact the technical support. |
|--|---|---|
| No. of pixels | 2 million pixels | 2 million pixels |
| Image transfer time | Ca. 34ms | Ca. 67ms |
| Dimensions (width x height x depth) | 29 x 29 x 58mm | 29 x 29 x 30mm |
| CCD size | 1/1.8 inch | 1/1.8 inch |
| Pixel size | 4.4 x 4.4µm | 4.4 x 4.4µm |

| | Quad-speed digital camera (ANPVC 1040D) | 5-Mega-Pixel digital Camera (ANPVC 1510D) |
|--|--|--|
| No. of pixels | 300,000 pixels | 5 million pixels |
| Image transfer time | Ca. 9ms | Ca. 66ms |
| Dimensions (width x height x depth) | 29 x 29 x 58mm | 44 x 44 x 57.5mm |
| CCD size | 1/3 inch | 2/3 inch |
| Pixel size | 7.4 x 7.4µm | 3.45 x 3.45µm |

| | Quad-speed Color digital camera (ANPVC 2040) | 2-Mega-Pixel Color digital Camera (ANPVC 2260) |
|--|---|---|
| No. of pixels | 300,000 pixels | 2 million pixels |
| Image transfer time | Ca. 9ms | Ca. 34ms |
| Dimensions (width x height x depth) | 29 x 29 x 58mm | 29 x 29 x 58mm |
| CCD size | 1/3 inch | 1/1.8 inch |

| | Quad-speed Color digital camera (ANPVC 2040) | 2-Mega-Pixel Color digital Camera (ANPVC 2260) |
|-------------------|---|---|
| Pixel size | 7.4 x 7.4 μ m | 4.4 x 4.4 μ m |

| | | |
|--|---|--|
| | Triple-speed Gray digital camera (ANPVC 5030) ➔ For details specification contact the technical support. | |
| No. of pixels | 300,000 pixels | |
| Image transfer time | Ca. 11ms | |
| Dimensions (width x height x depth) | 29 x 29 x 58mm | |
| CCD size | 1/3 inch | |
| Pixel size | 6.0 x 6.0µm | |

1.1.2 Connecting Cameras to the Imagechecker


To hook up a camera, connect it to the frame grabber camera connector via the camera cable (PoCL). Please read the notes below prior to operating any of the cameras.



◆ NOTE

- **Only use the cameras and camera cables recommended by Q.VITEC.**
- **Please refer to our homepage for a lens selection table:**
<http://www.qvitec.de>
- **Do not twist or bend the camera cables forcibly or apply load to the connector joints.**
- **Do not touch the CCD element or the lens surface of the camera. Attach the lens cap to keep off dust when not in use.**
- **When removing the cable from the connector, be sure to hold it by its plug to prevent excess force applying to the cable.**
- **Do not touch the terminals inside the connector of cameras and camera cables and take care not to allow foreign objects to come into the connectors.**

1.2 2-Mega-Pixel Digital Camera ANPVC 1210 / 2260



| No. | Description | Additional information |
|-----|-----------------|---|
| 1 | Lens | Use the proper lens for each camera. Make sure to use only lenses with a resolution suitable for the camera. For more information please contact Q.VITEC. The following lenses are available: <ul style="list-style-type: none"> • ANPVL162, lens f=16 for 2-mega-pixel digital camera • ANPVL252, lens f=25 for 2-mega-pixel digital camera • ANPVL502, lens f=50 for 2-mega-pixel digital camera |
| 2 | Lens fixture | C-Mount lens |
| 3 | Mounting plate | Insulation type plate. The plate can be attached to four sides of the camera. |
| 4 | Cable connector | Connects the camera cable Mini CL SDR (PoCL). |



◆ NOTE

- From now on, the 2-mega-pixel digital camera ANPVC1210 will be referred to as "2-mega camera".
- All camera settings like shutter speed, offset, gain, gamma, flash signal output etc. are made on the property page "Camera Settings" under Camera → Properties in the software.


1.2.1 Time Diagram for Image Grabbing

The camera ANPVC 1210 is a black-and-white progressive scan CCD camera. The camera is suitable for grabbing images of moving objects. The image grabbing is triggered by a signal from the frame grabber. With a delay of only a few μs the sensor will be completely exposed according to the shutter speed (exposure time) you have set. Immediately afterwards the frame grabber starts the transfer of the digital image signal to the memory of the Imagechecker, where it will be processed and evaluated.

There are 30 different settings between 1/15 seconds and 1/12500 seconds available for the shutter speed. The transfer time for an image from the camera to the memory is about 34ms for a full image with 1200 lines. This means a maximum image rate of about 30 full images per second.

When you are using partial scan mode, only a reduced number of image lines will be transferred and this reduces the transfer time. This helps you to increase the image or inspection

rate considerably. As the transfer time is not reduced proportionally to the number of image lines used, you need to determine the transfer time case by case.




Time diagram for image grabbing



◆ NOTE

Important: The duration of the image grabbing depends to a large degree on the image transfer time and the shutter speed:
image grabbing = shutter speed + image transfer time

All camera settings like shutter speed, offset, gain, gamma, flash signal output etc. are made on the property page "Camera Settings" under **Camera** → **Properties** in the software. For example, it is possible to set a delay between the camera trigger and the flash signal output so that you can set the optimum exposure time and duration for the image grabbing. This will help to reduce the interfering effect of other light sources and increase the service life of your flash. Make sure to activate the output of the flash signal with the option "Send flash signal".




Example for flash signal and camera trigger



◆ NOTE

Please also refer to the information supplied by the manufacturer of your flash!


1.2.2 Dimensions




1.2.3 Specifications

| Feature | | Specification |
|----------------------------|----------------------------|---|
| Image sensor | | Interline transfer method. 1/1.8" CCD, pixel size: 4.4μm x 4.4μm (tetragonal pixel) |
| No. of effective pixels | | 1628 x 1236 pixels (horizontal x vertical) |
| Grabbing mode | | Progressive scan, full image mode/partial scan mode |
| Shutter speed | | 1/15 - 1/12500s |
| Synchronous method | | Internal synchronous |
| Image output | | Digital output, adapted for progressive mount, 30fps |
| Lens fixture | | C-Mount lens |
| Rating | Power supply voltage | 12V DC +/- 10% |
| | Power consumption | 360mA maximum |
| Environmental requirements | Operational temperature | 0°C to +40°C (avoid ice and condensation) |
| | Operating ambient humidity | 35% to 85% RH (avoid ice and condensation at 25°C) |
| | Vibration resistance | 10 to 55 Hz, 1 sweep per minute, amplitude 1mm, 30 minutes each in X/Y/Z direction |
| | Shock resistance | 700 m/s ² , 3 times each in X/Y/Z direction |
| Weight | | About 65 g (without lens) |

1.2.4 Spectral Response



1.3 Quad-Speed Digital Camera ANPVC 1040



| No. | Description | Additional information |
|-----|-----------------|---|
| 1 | Lens | Use the proper lens for each camera. |
| 2 | Lens fixture | C-Mount lens |
| 3 | Mounting plate | Insulation type plate. The plate can be attached to four sides of the camera. |
| 4 | Cable connector | Connects the camera cable Mini CL SDR (PoCL). |



◆ NOTE


- **From now on, the quad-speed digital camera ANPVC1210 will be referred to as "quad-speed camera".**
- **All camera settings like shutter speed, offset, gain, gamma, flash signal output etc. are made on the property page "Camera Settings" under Camera → Properties in the software.**

1.3.1 Time Diagram for Image Grabbing

The camera ANPVC 1040 is a black-and-white progressive scan CCD camera. The camera is suitable for grabbing images of moving objects. The image grabbing is triggered by a signal from the frame grabber. With a delay of only a few μs the sensor will be completely exposed according to the shutter speed (exposure time) you have set. Immediately afterwards the frame grabber starts the transfer of the digital image signal to the memory of the Imagechecker, where it will be processed and evaluated.

There are 30 different settings between 1/10 seconds and 1/12500 seconds available for the shutter speed. The transfer time for an image from the camera to the memory is about 9ms for a full image with 480 lines. This means a maximum image rate of about 120 full images per second.

When you are using partial scan mode, only a reduced number of image lines will be transferred and this reduces the transfer time. This helps you to increase the image or inspection rate considerably. For example, if you only use two image lines and set a shutter speed of 100 μs , image grabbing only takes 2ms. This enables an image rate of 500 images per second. As the transfer time is not reduced proportionally to the number of image lines used, you need to determine the transfer time case by case.




Time diagram for image grabbing



◆ NOTE

Important: The duration of the image grabbing depends to a large degree on the image transfer time and the shutter speed:
image grabbing = shutter speed + image transfer time

All camera settings like shutter speed, offset, gain, gamma, flash signal output etc. are made on the property page "Camera Settings" under **Camera** → **Properties** in the software. For example, it is possible to set a delay between the camera trigger and the flash signal output so that you can set the optimum exposure time and duration for the image grabbing. This will help to reduce the interfering effect of other light sources and increase the service life of your flash. Make sure to activate the output of the flash signal with the option "Send flash signal".




Example for flash signal and camera trigger



◆ NOTE

Please also refer to the information supplied by the manufacturer of your flash!


1.3.2 Dimensions




1.3.3 Specifications

| Feature | | Specification |
|----------------------------|----------------------------|--|
| Image sensor | | Interline transfer method. 1/3" CCD, pixel size: 7.4μm x 7.4μm (tetragonal pixel) |
| No. of effective pixels | | 659 x 494 pixels (horizontal x vertical) |
| Grabbing mode | | Progressive scan, full image mode/partial scan mode |
| Shutter speed | | 1/10s - 1/12500s |
| Lens fixture | | C-Mount lens |
| Synchronous method | | Internal synchronous |
| Image output | | Digital output, adapted for progressive mount, 120fps |
| Rating | Power supply voltage | 12V DC +/- 10% |
| | Power consumption | 235mA maximum |
| Environmental requirements | Operational temperature | 0°C to +45°C (avoid ice and condensation) |
| | Operating ambient humidity | 35% to 85% RH (avoid ice and condensation at 25°C) |
| | Vibration resistance | 10 to 55 Hz, 1 sweep per minute, amplitude 1mm, 30 minutes each in X/Y/Z direction |
| | Shock resistance | 700 m/s ² , 3 times each in X/Y/Z direction |
| Weight | | About 65g (without lens) |

1.3.4 Spectral Response



1.4 5-Mega-Pixel Digital Camera ANPVC 1510



| No. | Description | Additional information |
|-----|-----------------------------|--|
| 1 | DC IN connector | For connecting the camera cable that also supplies power to the camera |
| 2 | Digital interface connector | For connecting the Camera Link cable |
| 3 | Camera cable | Mini CL-SDR (PoCL) |
| 4 | Camera Link cable | |
| 5 | Fastening screws | |



◆ NOTE


- **From now on, the 5-mega-pixel digital camera ANPVC 1510 will be referred to as "5-mega camera".**
- **All camera settings like shutter speed, offset, gain, gamma, flash signal output etc. are made on the property page "Camera Settings" under Camera → Properties in the software.**

1.4.1 Time Diagram for Image Grabbing

The camera ANPVC 1510D is a black-and-white progressive scan CCD camera. The camera is suitable for grabbing images of moving objects. The image grabbing is triggered by a signal from the frame grabber. With a delay of only a few μs the sensor will be completely exposed according to the shutter speed (exposure time) you have set. Immediately afterwards the frame grabber starts the transfer of the digital image signal to the memory of the Imagechecker, where it will be processed and evaluated.

There are 30 different settings between 1/10 seconds and 1/10000 seconds available for the shutter speed. The transfer time for an image from the camera to the memory is about 66ms for a full image with 2050 lines. This means a maximum image rate of about 15 full images per second.

If you would like to use the partial scan mode, the camera needs to be programmed via the serial interface. Please contact the technical support.




Time diagram for image grabbing



◆ NOTE

Important: The duration of the image grabbing depends to a large degree on the image transfer time and the shutter speed:
image grabbing = shutter speed + image transfer time

All camera settings like shutter speed, offset, gain, gamma, flash signal output etc. are made on the property page "Camera Settings" under **Camera** → **Properties** in the software. For example, it is possible to set a delay between the camera trigger and the flash signal output so that you can set the optimum exposure time and duration for the image grabbing. This will help to reduce the interfering effect of other light sources and increase the service life of your flash. Make sure to activate the output of the flash signal with the option "Send flash signal".




Example for flash signal and camera trigger



◆ NOTE

Please also refer to the information supplied by the manufacturer of your flash!

1.4.2 Dimensions




Unit: mm (inch)

1.4.3 Specifications

| Feature | | Specification |
|----------------------------|----------------------------|---|
| Image sensor | | Interline transfer method. 2/3" CCD, pixel size: 3,45µm x 3.45µm (tetragonal pixel) |
| No. of effective pixels | | 2456 x 2058 pixels (horizontal x vertical) |
| Grabbing mode | | Progressive scan, full image mode/partial scan mode |
| Shutter speed | | 1/10s – 1/10000s |
| Synchronous method | | Internal synchronous |
| Image output | | Digital output, adapted for progressive mount, 66fps |
| Lens fixture | | C-Mount lens |
| Rating | Power supply voltage | 12V DC +/- 10% |
| | Power consumption | maximal 3.6W |
| Environmental requirements | Operational temperature | -5 to +45°C (avoid ice and condensation) |
| | Operating ambient humidity | 20 to 80% relative humidity (avoid ice and condensation at 25°C) |
| | Vibration resistance | 20 to 55Hz, 1 sweep per minute, amplitude 1mm, 30 minutes each in X/Y/Z direction |
| | Shock resistance | 700m/s ² , 3 times each in X/Y/Z direction |
| Weight | | Ca. 135g (without lens) |


1.4.4 Spectral Response



1.5 Partial Scan Mode

Here you can specify which part of the image sensor should be used for image grabbing. If you reduce the sensor area used, the image grabbing time will be shorter. The image size for a full image is 640 x 480 pixels for the quad-speed and the ultra-compact camera, and 1600 x 1200 pixels for the 2-mega camera.

The image part to be captured is defined by a start line and an end line. The example shows a partial image to be captured from start point at 160 and end point at 320.



◆ NOTE

- **You can set up to two partial areas per camera.**
- **The following minimum values apply when you set the capturing area:**
1 line for the quad-speed camera ANPVC 1040 and the ultra compact camera ANPVC 1021.
100 lines for the 2 Mega camera ANPVC 1210.

Index

5

5-mega-pixel camera ANPVC 1510 20

A

ANPC 1510 20

 Dimensions 22

 Specifications 23

 Spectral response 23

 Time response 20

ANPVC 1040 16

 Dimensions 14

 Specifications 18

 Spectral response 19

 Time response 16

ANPVC 1210 12

 Dimensions 14

 Specifications 14

 Spectral response 15

 Time response 12

C

Camera specifications 14, 18, 23

Camera types 8

Connect camera 11

P

P400MD

 Connect camera 11

Partial scan mode 24

Q

Quad-speed digital camera ANPVC 1040 16

S

Specifications

 5-mega-pixel digital camera 23

 Megapixel camera 14

Spectral response 15, 19, 23

T

Time response 12, 16, 20

Record of Changes

| Manual No. | Date | Description of changes |
|-------------------|------------|------------------------|
| MAN-HW-CLCAM-15CE | March 2015 | First edition |